

## Experiment 2

# The Michelson Interferometer

### 2.1 Outline of Theory

In the Michelson Interferometer (Figure 2.1), light from an extended source  $S$  is split into two beams of equal intensity by the lightly silvered back surface of the beam splitter  $G_1$ .





Figure 2.1: The Michelson Interferometer – Theoretical Representation

One beam is reflected from the movable mirror  $M_1$ , whose position is determined by the micrometer screw  $V$ . The other beam passes through the compensator plate  $G_2$  (equal in optical thickness to  $G_1$ ), and is reflected from the adjustable mirror  $M_2$ . This beam is again reflected by  $G_1$ , and the interference phenomena between the two beams can be observed in the direction  $O$ .

### Fringes of Equal Inclination

When  $M_2$  is adjusted perpendicular to  $M_1$ , then its image  $M_2'$ , formed by reflection in  $G_1$  is parallel to  $M_1$  and separated by a distance  $d = M_2G_1 - M_1G_1$ . The source plane  $S$  gives, by reflection in  $G_1$ ,  $M_1$  and  $M_2$ , two virtual source planes  $S_1$  and  $S_2$  separated by  $2d$ .

The path difference  $\Delta$  between two rays, coming from  $P$  on  $S$ , at an angle  $\psi$  with the axis perpendicular  $M_1$ , is then



$$\Delta = 2d \cos \psi$$

The phase difference between these two rays is then

$$\delta = \left( \frac{2\pi}{\lambda} \right) 2d \cos \psi - \pi$$


since the returning beam from  $M_2$  suffers a phase change  $\pi$  on reflection from  $G_1$ . Looking in the direction  $G_1M_1$ , an observer will therefore see a set of circular fringes. The spacing of these fringes of equal inclination will depend on the value of  $d$ . The central portion of the field will be dark for

$$d = \frac{m\lambda}{2}, \quad m = 0, \pm 1, \pm 2, \dots$$

These fringes remain visible for relatively large values of  $d$ , depending on the monochromaticity of the source  $S$ .

### Fringes of Equal Thickness

When  $M_2$  is not exactly perpendicular to  $M_1$ , and  $d$  is small, then the space between  $M_1$  and  $M_2'$  is a thin wedge. The two rays from  $P$  now seem to diverge from a point  $P'$  near  $M_1$ .



For various positions of  $P$  on  $S$ , the point  $P'$  stays always close to  $M_1$  and the path difference between the two rays is constant. Looking at  $M_1$  we will see a set of straight fringes, parallel to the wedge edge. These fringes of equal thickness connect points of equal air thickness in the wedge formed by  $M_1$  and  $M_2'$ . If  $d$  increases, the fringes will become curved, and for large  $d$  no fringes can be seen.

### White Light Fringes

If  $S$  is a white light source, then fringes can only be seen for very small values of  $d$ . For  $d = 0$  the central fringe is black, the neighbouring fringes are strongly coloured, and about 5-10 fringes can be seen. If  $d$  is larger than a few wavelengths, then no fringes at all can be seen with a white light source.

## 2.2 Procedure

### 2.2.1 Alignment of $M_2$

1. Illuminate the diffusing screen  $S$  with the mercury lamp, using the green filter to obtain monochromatic light with  $\lambda = 5461 \text{ \AA}$ .
2. Place the small pointer on the screen  $S$ . Looking into the interferometer from  $O$ , you will see two images of this pointer, formed by reflection in  $M_1$  and  $M_2$ .
3. Adjust  $M_2$  until these images coincide. Circular fringes will then be seen behind  $M_1$ .

### 2.2.2 Calibration of the screw $V$ in terms of $\lambda$

Rotating the micrometer screw  $V$  displaces  $M_1$  in the direction of  $G_1M_1$ . Displacing  $M_1$  over a distance  $\lambda/2$  brings a new fringe into the field of view (the order of interference changes from  $m$  to  $m + 1$ ).

1. Measure the rotation of  $V$  necessary to change  $m$  by fifty, and repeat this measurement five times.

You can now calculate how far one full revolution of  $V$  displaces  $M_1$ , i.e. you have calibrated the effective pitch of  $V$  in terms of  $\lambda$ .

### 2.2.3 Measurement of the refractive index

1. Now that  $V$  is calibrated, make the distance  $G_1M_1$  closely equal to  $G_1M_2$ , by rotating  $V$  until you observe only one or two circular fringes.
2. Adjust  $M_2$  until fringes of equal thickness are visible.
3. Illuminate  $S$  with white light and very carefully adjust  $V$  until coloured (straight) fringes are observed. If the black fringe is in the centre of the field of view, then  $d$  is exactly zero. Let the reading of  $V$  in this position be  $q_1$ .
4. Place a microscope cover glass carefully against  $M_2$  so that it covers the lower half of  $M_2$ . The optical path length of that part of the beam going through the cover glass to and from  $M_2$  has now increased by  $2(n - 1)\ell$ , where  $n$  is the refractive index and  $\ell$  the thickness of the cover glass.
5. Increase the optical path  $G_1M_1G_1$  by an equal amount by displacing  $M_1$  with the help of  $V$  so that the coloured fringes will be seen in the lower half of the field of view.
6. Let the reading of  $V$  be  $q_2$  when the black fringe, as seen through the cover glass, is in the center of the field of view.

7. From  $q_1$  and  $q_2$  (each measured five times), and the calibration of  $V$  obtained in section 2.2.2, you can calculate  $n\ell$ .
8. Measure  $\ell$  five times with a micrometer, and calculate  $n$  (assuming  $n_{\text{air}}$  equal to 1.000).
9. Calculate the experimental error in your value of  $n$ .




Figure 2.2: The Michelson Interferometer – Physical Representation